

Features

- Operates down to 1.4 K with appropriate sensor
- One sensor input
- Supports diode and RTD sensors
- 0 V to 10 V or 4 mA to 20 mA output
- Large 5-digit LED display
- RS-232C serial interface and alarm relays

Model 211 Temperature Monitor



Product Description

The Lake Shore single-channel Model 211 Temperature Monitor provides the accuracy, resolution, and interface features of a benchtop temperature monitor in an easy to use, easily integrated, compact instrument. With appropriate sensors, the Model 211 measures temperature from 1.4 K to 800 K including temperatures in high vacuum and magnetic fields. Alarms, relays, user-configurable analog voltage or current output, and a serial interface are standard features on the Model 211. It is a good choice for liquefied gas storage and monitoring, cryopump control, cryo-cooler, and materials science applications, and for applications that require greater accuracy than thermocouples allow.

Sensor Input Reading Capability

The Model 211 Temperature Monitor supports diode temperature sensors and resistance temperature detectors (RTDs). The Model 211 can be configured for the type of sensor in use from the instrument front panel. Ensuring high accuracy and 5-digit measurement resolution are 4-lead differential measurement and 24-bit analog-to-digital conversion.


The Model 211 converts voltage or resistance to temperature units based on temperature response curve data for the sensor in use. Standard temperature response curves for silicon diodes and platinum RTDs are included in instrument firmware. The Model 211 also provides non-volatile memory for one 200-point temperature response curve, which can be entered via the serial interface.

Interface

With an RS-232C serial interface and other interface features, the Model 211 is valuable as a stand-alone monitor and is easily integrated into other systems. Setup and every instrument function can be performed via serial interface or the front panel of the Model 211. Temperature data can be read up to seven times per second over computer interface; the display is updated twice each second. High and low alarms can be used in latching mode for error limit detection and in non-latching mode in conjunction with relays to perform simple on-off control functions. The analog output can be configured for either 0 to 10 V or 4 to 20 mA output.

Display

The Model 211 has a 6-digit LED display with measurements available in temperature units K, °C, °F, or sensor units V or Ω.



- ① Power input connector
- ② Serial (RS-232C) I/O (DTE)
- ③ Analog output

Sensor Selection

Sensor Temperature Range (sensors sold separately)

		Model	Useful Range	Magnetic Field Use
Diodes	Silicon Diode	DT-670-SD	1.4 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-670E-BR	30 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-414	1.4 K to 375 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-421	1.4 K to 325 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-470-SD	1.4 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-471-SD	10 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	GaAlAs Diode	TG-120-P	1.4 K to 325 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
	GaAlAs Diode	TG-120-PL	1.4 K to 325 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
	GaAlAs Diode	TG-120-SD	1.4 K to 500 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
Positive Temperature Coefficient RTDs	100 Ω Platinum	PT-102/3	14 K to 873 K	$T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$
	100 Ω Platinum	PT-111	14 K to 673 K	$T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$
	Rhodium-Iron	RF-800-4	1.4 K to 500 K	$T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$
	Rhodium-Iron	RF-100T/U	1.4 K to 325 K	$T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$
Negative Temperature Coefficient RTDs¹	Cernox™	CX-1010	2 K to 325 K ⁴	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1030-HT	3.5 K to 420 K ^{2,5}	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1050-HT	4 K to 420 K ^{2,5}	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1070-HT	15 K to 420 K ²	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1080-HT	50 K to 420 K ²	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Germanium	GR-200A/B-1000	2.2 K to 100 K ³	Not Recommended
	Germanium	GR-200A/B-1500	2.6 K to 100 K ³	Not Recommended
	Germanium	GR-200A/B-2500	3.1 K to 100 K ³	Not Recommended
	Carbon-Glass	CGR-1-500	4 K to 325 K ⁴	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Carbon-Glass	CGR-1-1000	5 K to 325 K ⁴	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Carbon-Glass	CGR-1-2000	6 K to 325 K ⁴	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Rox™	RX-102A	1.4 K to 40 K ⁴	$T > 2 \text{ K} \ \& \ B \leq 10 \text{ T}$

Silicon diodes are the best choice for general cryogenic use from 1.4 K to above room temperature. Diodes are economical to use because they follow a standard curve and are interchangeable in many applications. They are not suitable for use in ionizing radiation or magnetic fields.

Cernox™ thin-film RTDs offer high sensitivity and low magnetic field-induced errors over the 2 K to 420 K temperature range. Cernox sensors require calibration.

Platinum RTDs offer high uniform sensitivity from 30 K to over 800 K. With excellent reproducibility, they are useful as thermometry standards. They follow a standard curve above 70 K and are interchangeable in many applications.

¹ Single excitation current may limit the low temperature range of NTC resistors

² Non-HT version maximum temperature: 325 K

³ Low temperature limited by input resistance range

⁴ Low temperature specified with self-heating error: ≤ 5 mK

⁵ Low temperature specified with self-heating error: ≤ 12 mK

Typical Sensor Performance – see Appendix F for sample calculations of typical sensor performance

	Example Lake Shore Sensor	Temp	Nominal Resistance/Voltage	Typical Sensor Sensitivity ⁶	Measurement Resolution: Temperature Equivalents	Electronic Accuracy: Temperature Equivalents	Temperature Accuracy including Electronic Accuracy, CalCurve™, and Calibrated Sensor
Silicon Diode	DT-670-SD with 1.4H calibration	1.4 K	1.644 V	-12.49 mV/K	1.6 mK	±26 mK	±38 mK
		77 K	1.028 V	-1.73 mV/K	11.6 mK	±152 mK	±174 mK
		300 K	0.5597 V	-2.3 mV/K	8.7 mK	±94 mK	±126 mK
		500 K	0.0907 V	-2.12 mV/K	9.4 mK	±80 mK	±130 mK
Silicon Diode	DT-470-SD-13 with 1.4H calibration	1.4 K	1.6981 V	-13.1 mV/K	1.5 mK	±26 mK	±38 mK
		77 K	1.0203 V	-1.92 mV/K	10.5 mK	±137 mK	±159 mK
		300 K	0.5189 V	-2.4 mV/K	8.4 mK	±88 mK	±120 mK
GaAlAs Diode	TG-120-SD with 1.4H calibration	1.4 K	5.391 V	-97.5 mV/K	0.2 mK	±13 mK	±25 mK
		77 K	1.422 V	-1.24 mV/K	16.2 mK	±359 mK	±381 mK
		300 K	0.8978 V	-2.85 mV/K	7 mK	±120 mK	±152 mK
100 Ω Platinum RTD 500 Ω Full Scale	PT-103 with 1.4J calibration	30 K	3.66 Ω	0.19 Ω/K	10.5 mK	±25 mK	±35 mK
		77 K	20.38 Ω	0.42 Ω/K	4.8 mK	±20 mK	±32 mK
		300 K	110.35 Ω	0.39 Ω/K	5.2 mK	±68 mK	±91 mK
		500 K	185.668 Ω	0.378 Ω/K	5.3 mK	±109 mK	±155 mK
Cernox™	CX-1050-SD-HT ⁷ with 4M calibration	4.2 K	3507.2 Ω	-1120.8 Ω/K	45 μK	±1.4 mK	±6.4 mK
		77 K	205.67 Ω	-2.4116 Ω/K	20.8 mK	±75.6 mK	±91.6 mK
		300 K	59.467 Ω	-0.1727 Ω/K	290 mK	±717 mK	±757 mK
Germanium	GR-200A-1000 with 1.4D calibration	420 K	45.03 Ω	-0.0829 Ω/K	604 mK	±1.43 K	±1.5 K
		2 K	6674 Ω	-9930 Ω/K	5 μK	±0.3 mK	±4.3 mK
		4.2 K	1054 Ω	-526 Ω/K	95 μK	±10 mK	±14 mK
		10 K	170.9 Ω	-38.4 Ω/K	1.3 mK	±4.4 mK	±9.4 mK
Carbon-Glass	CGR-1-2000 with 4L calibration	100 K	2.257 Ω	-0.018 Ω/K	2.78 K	±5.61 K	±5.77 K
		4.2 K	2260 Ω	-2060 Ω/K	25 μK	±0.5 mK	±4.5 mK
		77 K	21.65 Ω	-0.157 Ω/K	319 mK	±692 mK	±717 mK
		300 K	11.99 Ω	-0.015 Ω/K	3.33 K	±7 K	±7.1 K

⁶ Typical sensor sensitivities were taken from representative calibrations for the sensor listed

⁷ Non-HT version maximum temperature: 325 K

Specifications

Input Specifications

	Sensor Temperature Coefficient	Input Range	Excitation Current	Display Resolution	Measurement Resolution	Electronic Accuracy
Diode	negative	0 V to 2.5 V	10 μA ±0.05% ⁸	100 μV	20 μV	±160 μV ±0.01% of rdg
	negative	0 V to 7.5 V	10 μA ±0.05% ⁸	100 μV	20 μV	±160 μV ±0.02% of rdg
PTC RTD	positive	0 Ω to 250 Ω	1 mA ±0.3% ⁹	10 mΩ	2 mΩ	±0.004 Ω ±0.02% of rdg
	positive	0 Ω to 500 Ω	1 mA ±0.3% ⁹	10 mΩ	2 mΩ	±0.004 Ω ±0.02% of rdg
	positive	0 Ω to 5000 Ω	1 mA ±0.3% ⁹	100 mΩ	20 mΩ	±0.06 Ω ±0.04% of rdg
NTC RTD	negative	0 Ω to 7500 Ω	10 μA ±0.05% ⁸	100 mΩ	50 mΩ	±0.1 Ω ±0.04% of rdg

⁸ Current source error has negligible effect on measurement accuracy

⁹ Current source error is removed during calibration

Thermometry

Number of inputs 1

Input configuration Input can be configured from the front panel to accept any of the supported input types

Isolation Measurement is not isolated from chassis ground

A/D resolution 24-bit

Input accuracy Sensor dependent – refer to Input Specifications table

Measurement resolution Sensor dependent – refer to Input Specifications table

Maximum update rate 7 rdg/s

User curve One 200-point CalCurve™ or user curve in non-volatile memory

Sensor Input Configuration

	Diode/RTD
Measurement type	4-lead differential
Excitation	Constant current
Supported sensors	Diodes: Silicon, GaAlAs RTDs: 100 Ω Platinum, 1000 Ω Platinum, Carbon-Glass, Cernox™, and Rox™
Standard curves	DT-470, DT-670, CTI-C, PT-100, and PT-1000
Input connector	Shared 25-pin D-sub

Front Panel

Display 5-digit LED
Number of reading displays 1
Display units K, °C, °F, V, and Ω
Reading source Temperature and sensor units
Display update rate 2 rdg/s
Temp display resolution 0.001° from 0° to 99.999°, 0.01° from 100° to 999.99°, 0.1° above 1000°

Sensor units

display resolution Sensor dependent to 5 digits
Display annunciators K, °C, °F, and V/Ω
Keypad 4 full travel keys, numeric and specific functions
Front panel features Display brightness control, keypad lock-out

Interface

Serial interface

Electrical format RS-232C
Max baud rate 9600 baud
Connector 9-pin D-sub
Reading rate Up to 7 rdg/s

Alarms

Number 2, high and low
Data source Temperature
Settings High setpoint, Low setpoint, Dead band, Latching or Non-latching
Actuators Display message, relays

Relays

Number 2
Contacts Normally Open (NO), Normally Closed (NC), and Common (C)
Contact rating 30 VDC at 1 A
Operation Activate relays on high or low input alarm or manual
Connector Shared 25-pin D-sub

Analog output

Isolation Output is not isolated from chassis ground
Update rate 7 readings per s
Data source Temperature

	Voltage	Current
Range	0 V to 10 V	4 mA to 20 mA
Accuracy	± 1.25 mV	± 2.5 μA
Resolution	0.3 mV	0.6 μA
Min load resistance	500 Ω	NA
Compliance voltage	NA	10 V
Load regulation	NA	± 0.02% of reading 0 to 500 Ω

Scales:	Temperature	Sensor units (fixed by type)
	0 K to 20 K	Diodes: 1 V = 1 V
	0 K to 100 K	100 Ω Platinum: 1 V = 100 Ω
	0 K to 200 K	1000 Ω Platinum: 1 V = 1000 Ω
	0 K to 325 K	NTC Resistor: 1 V = 1000 Ω
	0 K to 475 K	
	0 K to 1000 K	

Settings Voltage or current, scale
Connector Shared 25-pin D-sub

General

Ambient temperature

Range 15 °C to 35 °C at rated accuracy,
10 °C to 40 °C at reduced accuracy

Power requirements

Size Regulated +5 VDC at 400 mA
96 mm W × 48 mm H × 166 mm D
(3.8 in × 1.9 in × 6.5 in)

Mounting

Panel mount into 91 mm W × 44 mm H
(3.6 in × 1.7 in) cutout

Weight

0.45 kg (1 lb)
CE mark, RoHS compliant

Power Supply (109-132)

Power requirements 100 – 240 VAC, 50 or 60 Hz, 0.3 A max
Output +5 V at 1.2 A
Size 40.5 mm W × 30 mm H × 64 mm D
(1.6 in × 1.2 in × 2.5 in)
Weight 0.15 kg (0.33 lb)



2111 Single 1/4 DIN panel-mount adapter, 105 mm W × 132 mm H (4.1 in × 5.2 in)



2112 Dual 1/4 DIN panel-mount adapter, 105 mm W × 132 mm H (4.1 in × 5.2 in)

Ordering Information

Part number	Description
211S	Model 211 temperature monitor, single channel
211N	Model 211S with no power supply

Accessories Included with 211S

109-132 100-240 V, 6 W power supply (universal input, interchangeable input plugs)
G-106-253 Sensor input mating connector (DB-25)
G-106-264 Shell for sensor input mating connector
 Calibration certificate
MAN-211 Model 211 user manual

Options and accessories

2115 Power supply splitter cable — allows two Model 211Ss to be powered from one supply
2111 Single 1/4 DIN panel-mount adapter
2112 Dual 1/4 DIN panel-mount adapter
8000 CalCurve™, CD-ROM (included with calibrated sensor)
8001-211 CalCurve™, factory installed
CAL-211-CERT Instrument recalibration with certificate
CAL-211-DATA Instrument recalibration with certificate and data

